

## QUOTATION

Number: 1433602 / 21/03/2023

To: ELTA 90 Medical Genetics LLC

In attn.:

Deer Partner,

With thanks for our interest to our products, we are sending you below our price offer:

Code	Product	UM	Quantity	Price / UM	Currency
MNA51000	<b>Stereomicroscope Nikon SMZ800N</b> Stereomicroscope Nikon SMZ800N - zoom ratio: 8:1 - zooming range: 1-8 - working distance: 78 mm - trinocular tube, with a port for C-mount Camera - interpupillary distance adjustment: 50 mm or wider - Eyepiece inclination: 10° - Eyepieces (with diopter adjustment): 10x(FN 22), with diopter adjuster, rubber eye-shields, reticle lead - objective: P-PLAN 1X, W.D. 78 mm	BUC	1.00	6,267.00	EUR

This configuration is ready to be installed in the workstation.

Equipped with photo-video documentation system Mshot MS60:

Effective pixels: 6.3 megapixels  
resolution: 3072\*2048  
Sensor size: 1/1.8 inches  
Pixel size: 2.4 um\*2.4 um  
Frame rate:  
71FPS □ 3072x2048  
89FPS: 1536×1024(BIN2)  
166FPS: 1920×1080(1080P)  
178FPS: 1280×720(720P)

Any size ROI  
Exposure time: 6us □ 7s  
Effective gain: 1 □ 16X  
Spectral response: 380nm □ 650nm  
Scanning mode: Progress scan / continuous output  
Trigger mode: Soft trigger  
Shutter: Electronic rolling shutter  
A / D conversion/: 12Bit  
Exposure: Manual/Auto/Area exposure  
White balance: Manual/One push/Area  
Image format: TIF,BMP,JPG,RAW  
Soft ware Interface: DirectShow / TWAIN  
Operating system: Windows XP 32Bit, Windows 7/8/10  
32/64Bit  
Data interface: USB3.0 5Gbps B type interface

**Total (without VAT) : 6,267.00 EUR**

**Obs:**

**Commercial terms:**

**Prices:** do not include VAT, but include all taxes and commissions on Romanian territory

**Validity of the offer:** the offer is valid until 19/06/2023

**Payment deadline:**

**Payment terms:**

**Account details:**

RO80BACX0000000533262003 - Unicredit Bank

Delivery time: up to 2 business days from the date of receipt of the order for the products in stock. For non-stock products, the delivery time will be communicated after your order confirmation.

To make it easier to send your order, you can use this as a form of order. You will receive our confirmation with the delivery details as soon as possible.

Hoping that our offer meets your expectations, we are waiting you order by email address:office@elta90mr.ro.

Drafted by the:

SC ELTA 90 MEDICAL RESEARCH SRL .

Agent:Corina Dragne

Phone agent:021.232.26.94/96

Email agent:service@elta90mr.ro

Client name: ELTA 90 Medical Genetics LLC

Endorsed by the

Date:

Signature and stamp: